

CALL US TODAY!
(610) 337-3762
For Promotional Pricing

IsoSciences offers the most extensive selection of Vitamin K internal standards on the market.

In addition to the d7 isotopically labeled versions, IsoSciences has developed the only $^{13}\text{C}_6$ labeled K2 standards available. If you are interested in a Vitamin K analyte not listed below, please contact us.

Vitamin K1

		CATALOG NUMBER	PURITY	ISOTOPE INCORPORATION
K1-[d7]	Phytonadione-[d7]	5072	≥97%	≥99%
K1-[$^{13}\text{C}_6$]	Phytonadione-[$^{13}\text{C}_6$]	14329	≥98%	≥99%

Vitamin K2

		CATALOG NUMBER	PURITY	ISOTOPE INCORPORATION
MK-4-[d7]	Menaquinone-4-[d7]	14331	≥95%	≥98%
MK-4-[$^{13}\text{C}_6$]	Menaquinone-4-[$^{13}\text{C}_6$]	15277	≥95%	≥99%
MK-7-[d7]	Menaquinone-7-[d7]	15103	≥95%	≥98%
MK-7-[$^{13}\text{C}_6$]	Menaquinone-7-[$^{13}\text{C}_6$]	15278	≥95%	≥99%
MK-9-[d7]	Menaquinone-9-[d7]	15406	≥95%	≥98%
MK-9-[$^{13}\text{C}_6$]	Menaquinone-9-[$^{13}\text{C}_6$]	15279	≥95%	≥99%

Vitamin K3

		CATALOG NUMBER	PURITY	ISOTOPE INCORPORATION
K3-[d8]	Menadione-[d8]	7229	≥97%	≥98%

Epoxides

		CATALOG NUMBER	PURITY	ISOTOPE INCORPORATION
K1 2,3-Epoxide	Phytonadione 2,3-Epoxide	16002	≥95%	N/A
K1-[d7] 2,3-Epoxide	Phytonadione-[d7] 2,3-Epoxide	16003	≥95%	≥98%
MK-4 2,3-Epoxide	Menaquinone-4 2,3-Epoxide	16004	≥95%	N/A
MK-4 [d7] 2,3-Epoxide	Menaquinone-4-[d7] 2,3-Epoxide	16005	≥95%	≥98%

MENAQUINONE-4-[$^{13}\text{C}_6$]

